

**Product name: Medi-Relief wax- Art.Nr. 200100****Swiss MediTec***Safety data sheet*

(REACH regulation (EC) n° 1907/2006 - n° 453/2010)

**Section 1. Identification of the Substance/Preparation and the Company/Undertaking**

Code: 200100

Trade Name: Orthocare Relief Wax

Supplier and Distributor: Swiss MediTec

Chemical Family: Microcrystalline Paraphin

Chemical Name: Formulated Product

Formula: N/A

**Section 2. Composition**

Hazardous Ingredients: None

CAS Number

Exposure Limits

Concentration

N/A

N/A

N/A

Hazardous materials identification system (HMIS)

(For material as packaged):

HMIS Ratings:

Health

Flammability

Reactivity

0

0

0

0= Minimal

1= Slight

2= Moderate

3= Serious

4= Severe

**Section 3. Hazards Identification****Hazardous Ingredients:**

This product has non-hazardous ingredient as defined under the criteria of the federal OSHA hazard communication standard 29 CFR 1910 1200

**Medical conditions****Aggravated by exposure:**

None

**Exposure limits:**

None

**Potential health effects:**

None

**Route(s) of Entry:**

Ingestion

Human effects and symptoms of overexposure:

None

**Acute inhalation:** N/A**Chronic inhalation:** None**Acute eye contact:** N/A**Chronic eye contact:** None**Carcinogenicity:** None**Section 4. First aid measures****First aid for eyes:**

In case of contact, flush eyes well with water for 15 minutes. Obtain medical attention if irritation occurs.

**First aid for skin:**

Remove excess material from the skin with a waterless skin cleaner. Flush skin with plenty of water and soap.

**First aid for inhalation:** N/A**First aid for ingestion:** Due to the physical state of this material, ingestion is unlikely to occur.**Section 5. Fire fighting measures****Flash point:** Not combustible**Flammable limits:** N/A**Upper Explosive Limit (UEL):** N/A**Lower Explosive Limit (LEL):** N/A**Extinguishing Media:** Foam, Dry chemicals, CO<sub>2</sub>, Water fog, Sand/Earth**Special fire fighting procedures:** N/A**Section 6. Accidental release measures****Spill or leak procedures:** Remove by mechanical means**Section 7. Handling and storage****Storage temperature (MIN/Max):** Ambient**Shelf life:** None**Special sensitive:** Do not store or mix with strong oxidants**Handling/storage precautions:** Store in sealed containers. Avoid excessive heat or direct sunlight.**Section 8. Exposure Controls****Eye protection requirements:** None

**Product name: Medi-Relief wax- Art.Nr. 200100****Swiss MediTec****Skin protection requirements:** None required**Ventilation requirements:** None**Respirator requirements:** None**Additional protective measures:** None**Section 9. Physical and chemical properties****Physical form:** Paste**Colour :**White ,Orange ,Yellow ,Purple, Green ,Red**Odour :**Odourless Orange Lemon Grape Lime Cherry**pH:** N/A**Boiling Point:** N/A**Melting-Freezing Point:** 140/150F (Melting)**Solubility in Water:** Insoluble**Specific Gravity:** 35/36G/CM3 (in Water=1)**Bulk Density:** N/A**% Volatile by Volume:** N/A**Vapour Pressure:** N/A**Section 10. Stability and reactivity****Stability:** Stable at ambient temperature and atmospheric pressure.**Hazardous Polymerisation:** Will not occur**Incompatibilities:** None Known**Decomposition Products:** By Fire; Carbon Monoxide, Carbon Dioxide, and Various Hydrocarbon Fragments.**Section 11. Toxicological information**

Toxicological testing has not been conducted with this material

**Section 12. Ecological information**

None

**Section 13 Disposal Considerations**

Incinerate

**Section 14. Transport information****Technical Shipping Name:** N/A**Freight Class Bulk:** N/A**Freight Class Package:** N/A**Product Label:** N/A**Hazard Class or Division:** Non-Hazardous**Hazard Class Division Number:** Not Hazardous by D.O.T. Regulations.**Section 15. Regulatory information****SARA Title III:****Section 302 Extremely****Hazardous Substances:** None**Section 311/312****Hazard Categories:** None**Section 313****Toxic Chemicals:** None**RERA Status:** Non-Hazardous**Section 16. Other Information**

This information is given in good faith, being based on the latest knowledge available to Swiss MediTec. No known relevant information has been omitted from the product information, risk information and safety precaution sheets. However, they shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship. The user should not assume therefore, the product is suitable for any abnormal use. If the information provided is insufficient to ensure safety in any particular application, please contact:

Swiss MediTec. prior to use for further advice